


TestWise Cup 5.20 Software


Data Collection from a Digital Force Gauge

Specially designed for applications where the whole package is peeled at once




Data collection window


Applications:



Food/Medical Cups



Thermoformed Trays and Blisters



Whole Pouches

Standard Features:

- Graphic representation of peel force vs. distance or force vs. time
- Calculates the average force between cursors, the initial and final peaks located outside the cursors and the area under the curve for each test
- Automatic detection of serial port and force gauge settings
- User selectable trigger method and length of data acquisition
- All data is saved in text format for compatibility with other software packages
- Works in both compression or tension modes
- User can change the position of the cursors for calculations
- Statistical calculations (\bar{x} , σ , $\bar{x}-2\sigma$, $\bar{x}-3\sigma$, min & max)
- Automatic field increment
- Supports printed report

Price List

TestWise Data Collection Software	\$ 875
Data Cable (required for data collection)	\$ 75
Overnight Shipping	\$ 40



Settings window



Vinatoru Enterprises, Inc.
209 West Hanover Road, Graham, NC 27253

tel. (336) 227-4300 fax (336) 227-5544
E-mail info@vinatoru.com
visit us at www.vinatoru.com