

# TestWise Pull 5.20 Software

## Data Collection from a Digital Force Gauge

Specially designed to peel a 1" wide sample according to ASTM F88



### Applications:



1" Wide Tray Samples



1" Wide Pouch Samples

Data collection window

### Standard Features:

- Graphic representation of peel force vs. distance or force vs. time
- Calculates the average force, the overall peak and the area under the curve for each test. The first and last 10% of the test are not taken into consideration
- The cursors are set at 10% and 90% of the test length
- Statistical calculations ( $\bar{x}$ ,  $\sigma$ ,  $\bar{x}-2\sigma$ ,  $\bar{x}-3\sigma$ , min & max)
- Automatic detection of serial port and gauge units of measurement
- Works in both compression or tension modes
- User selectable trigger method and length of the data acquisition
- All data is saved in text format for compatibility with other software packages
- Automatic field increment
- Supports printed report



Settings window

### Price List

TestWise Data Collection Software	\$ 875
Data Cable (required for data collection)	\$ 75
Overnight Shipping	\$ 40



Vinatoru Enterprises, Inc.  
209 West Hanover Road, Graham, NC 27253

tel. (336) 227-4300 fax (336) 227-5544  
E-mail info@vinatoru.com  
visit us at www.vinatoru.com